




Model Chart



Available sizes:

		WIDTH (WW)				
inch		16	17	18	19	20
cm		40	43	45	48	50
DEPTH (DD)	16					
	17					
	18					
	19					
	20					

width x depth

cm	inch
40 x 40	16 x 16
43 x 43	17 x 17
45 x 45	18 x 18
48 x 48	19 x 19
50 x 50	20 x 20

Label colour

	Blue
	Yellow
	Red
	Green
	Orange

Available Size Not Available

2 weeks lead time

Additional size information:

In the table below you can find the wheelchair seat depths and seat widths that the available cushion sizes can be used for. These dimensions are marked in the cushion's label colour.

		WHEELCHAIR SEAT WIDTH (WW)																	
inch		14	14 ¼	15	15 ½	16	16 ½	16 ¾	17	17 ½	18	18 ½	18 ¾	19	19 ½	20	20 ½	20 ¾	
cm		36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	
WHEELCHAIR SEAT DEPTH (DD)	15 ½																		
	16																		
	16 ½																		
	16 ¾																		
	17																		
	17 ½																		
	18																		
	18 ½																		
	18 ¾																		
	19																		
	19 ½																		
	20																		
	20 ½																		
20 ¾																			



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au