

rehabhire & sales

MEYRA<sup>®</sup> NOVIS

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

## XXL & XXL HD

The Chair for Rugged Types



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

## Tough and Stable

The Meyra XXL and XXL HD don't back down. Built to move bariatric users, with a great, stable and sturdy design, the XXL and XXL HD are tough. Thanks to the power tilt-in-space and recline function, as well as the choice of powered separate or centre mounted footrests, you can be tough and comfortable.

## Features and Benefits

- Ideal for those with larger body frames, supporting up to 200kg: XXL | 250kg: XXL HD
- Equipped with powerful motors for fast and efficient mobility: 2 x 300W or 2 x 350W motor options
- Power seat tilt and power backrest recline options available for pressure relief and health benefits
- Rear Wheel Drive for a smooth ride with larger casars and high top speed
- Built with reinforced components for strength and stability
- Wide range of standard and power options, including power tilt-in-space and recline function, and powered separate or centre mounted footrests
- LED lighting system for safe navigation in low light conditions
- Customisable with five frame and accent colours



Reinforced and built tough



Powerful 2 x 300W or  
2 x 350W motors



Powered tilt-in-space



Long-lasting  
LED lighting system

## Specifications

Length	1150mm with legrests   850mm without legrests
Width	670mm
Height without Headrest	1020mm - 1150mm
Transport Length	850mm
Transport Height Minimum	650mm
Seat Width	530mm - 700mm
Seat Depth	430mm - 600mm
Seat Height   Front with   without cushion	420mm - 470mm   460mm - 560mm   530mm - 630mm
Armrest Height	170mm - 280mm
Leg Length Adjustment	370mm - 510mm
Wheels   Front	260mm x 70mm
Wheels   Rear	356mm x 75mm
Mechanical Tilt in Space	0° to 10°
Electric Tilt in Space	0° to +18°
Mechanical Backrest Recline	-10° to +30°
Power Backrest Recline	-10° to +50°
Turning Radius	950mm
Permitted Gradient	13°
Maximum Obstacle Height	60mm
Range	30km - 35km 70Ah batteries
Driving Speeds	6kmph   10kmph
Motor Options	2 x 300W   2 x 350W
Empty Transport Weight	100kg
SWL	XXL: 200kg   XXL HD: 250kg





1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au