

reabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



iQ

Sets the new standard



1300 000 030 | reabhire.com.au

FWD Front Wheel Drive

Large wheels at front for better climbing of obstacles and tighter turning

- Front chassis clearance for easier transfer and drive wheels positioned at front to move in closer to tables, sinks, and cupboards.
- Short wheel base 82cm for easier use in confined space.



MWD Mid Wheel Drive

Highly manoeuvrable with the tightest turning circle in the range. Also features unique suspension

- Extra set of castors and central drive wheels for excellent stability and smoother drive comfort on uneven terrain.



RWD Rear Wheel Drive

Positioning of the wheel improves power and torque

- With only 4 wheels and the majority of the weight in the rear, this provides improved power and torque for difficult terrain.



Improved functionality for indoor use

Lowest floor to seat plate height 38cm with a standard seat elevation of 30cm allowing better access under tables, desks and when driving vehicles.

The TA-iQ MWD offers a 45cm small turning radius, ideally suited for indoor use to allow improved access in bathrooms and smaller spaces

Highest standard elevate 38cm to 68cm to allow for improved access to lower cupboards and higher benches

Improved performance for outdoor use

Offering a superior top speed of 10km/h and an improved range of 40km, with a softer dynamic suspension the TA iQ handles rougher terrain and climbs obstacles easily with improved rider comfort. This may also reduce pain and the risk of spasm for some users.

Superior safety when travelling by car

Crash tested to international safety standards ISO7176-19 and ISO 10542-5 whilst the occupant remains seated in the chair.

Specifications

Length without Legrests	900mm
Width without Seat	630mm
Minimum Height Floor to Seat Base	380mm
Minimum Turning Radius	FWD: 650mm MWD: 450mm RWD: 650mm
Electric Lift Height	300mm 380mm - 680mm
Motors	4-pole
Electric Tilt	45°
Maximum Range	40km
Maximum Speed	10km/h
SWL	140kg
Warranty	2 years (drive wheels excluded)

Notes



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au