

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

ELATA®

# Home Care Mattress



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

## Why we exist

A pressure care mattress is one intervention in a holistic plan. Always assess the whole patient/client.

Elata believe in zero preventable pressure injuries. It's ambitious but it's necessary. Pressure injuries cause complications. They also add workload and cost. We act before that pain starts.

## What this means

- Care that preserves dignity.
- Compassion that simplifies decisions.
- Commitment to better outcomes, every day.

## How we deliver

Advanced systems that adapt to every patient and pathway. Designed to simplify prescription, enhance comfort, and reduce clinical risk.

We deliver personalised, intuitive pressure care, protecting life with dignity, comfort, and trust.

## Every Pressure Injury Risk Level. One Matching Solution.

### Simple Choices. Clear Pathways.

#### Elata® Range. Quick Category Guide


Home Care	E100   E200	Simple prevention for home and aged care. Default choice in home care pathways.	N/A
Clinical Foam	E300	Non-powered. Prevents and treats Stage 1–2 pressure injuries. 9-zone foam with side wall support and heel offloading. Low maintenance.	
Clinical Hybrid	E400	Use alone for Stage 2 pressure injuries, or step up with a pump to support treatment up to Stage 4. Ideal for home and palliative care. Converts as needs change.	EP10AU
Clinical Alternating Hybrid	E500	Alternating air/foam system. Combines comfort with therapy. Treats up to Stage 4 pressure injuries.	EP20AU
Clinical Alternating Air	E600	Full alternating air (1-in-2 configuration). Treats up to Stage 4 pressure injuries. Ideal for surgical wards, acute wards, or clients with frailty/reduced mobility.	EP10AU
	E700	Advanced alternating (configurable 1-in-4). Auto-setup. Treats up to Stage 4 pressure injuries. For ICU, stable spinal, complex rehab, immobile clients.	EP20AU

The EP20AU Pump is fully backwards compatible. It can be used with all Elata pump-ready mattresses (E400 to E700).

## E100

### Essential Foam Care

Combines comfort foam, high density support foam, and a stable hard base. With side wall reinforcement, cooling fabric, and anti-slip backing, it delivers safe, reliable everyday care for low-risk clients.



### SWL

160kg

### Compliance

The E100 and E200 mattress fabric and foam meet the OEKO-TEX standard 100. The foam is also tested to be free from heavy metals, phthalates, CFC, and HCFC.


### Warranty

5 years

## E200

### Enhanced with Pocket Springs

Builds on the E100 with an integrated pocket spring core beneath its comfort foam. The upgrade delivers cooler airflow, adaptive contouring and greater stability for safer transfers. With the same cooling cover, high density foams, and strong side wall support as the E100, the E200 is the choice for advanced comfort and protection at home.



**S** **KS**

Cannot be used with a pump

### SWL

160kg

### Compliance

The E100 and E200 mattress fabric and foam meet the OEKO-TEX standard 100. The foam is also tested to be free from heavy metals, phthalates, CFC, and HCFC.

### Warranty

5 years





1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au