

rehabhire & sales

novis

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

## Netti 4U CE Plus - Tilt in Space and Reclining Wheelchair



*inspire joy of life*

**Netti**<sup>®</sup>  
4U CE Plus

1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

# Netti 4U CE Plus - Tilt in Space and Reclining Wheelchair

The Netti 4U CE Plus combines economy and individuality at the highest level. The frame design and extensive range of accessories meet and often exceed user demands. Well specified as standard and highly adaptable with a range of accessories, the Netti 4U CE Plus is a flexible and multi-functional wheelchair and offers the following features as standard:

## Flexible Back Solution

The Netti 4U CE Plus is equipped with a one-piece velcro back and a back cushion that can be adjusted for the user's height. Furthermore, you can choose from our complete product range of cushions, giving different levels of support for the trunk.

## The Narrowest Wheelchair on the Market

The overall width of the wheelchair is only 18 cm wider than the seat width with standard wheels.

## Folding Brakes

Folding brakes for easy and safe sideways and /or sling transfer.



## Hygienic Cover

This Netti 4U CE Plus is supplied with hygienic Easy Care covers which are wipeable. This makes the chair well suited for cases of incontinence as well as for multiple residents in hospitals and facilities.

## Seat Width Flexibility

Comfort pads and seat width extension options ensure an economic solution as well as an ideal fit for the user. The seat width can be adjusted 20 or 40 mm in or out, using the seat width extension kit.

## One-Hand-Operated Durable Arm Support

The functional designed arm support offers increased stability to the user. They are height and depth adjustable and easy to mount and dismount with one hand, making it great for sideways or sling transfer.




350-500 mm options

Seat width



425-500 mm

Seat depth



500 mm (600 mm)

Back height



465-500 mm

Seat height



-5° - +20°

Tilt



90-135°

Back angle



135 kg

Max load



27,5-29 kg

Weight



Max. 135 kg

Crash tested

4U CE Plus Tilt in Space Wheelchair	35 cm Seat Width	N-D64035
	40 cm Seat Width	N-D64040
	45 cm Seat Width	N-D64045
	50 cm Seat Width	N-D64050

## Available as a standard configuration suitable to most users' needs

Well specified standard features and configuration make the Netti 4U CE Plus a great option to suit most users' needs. This makes it a cost-effective seating solution without compromising on quality.

### Standard Specifications:

Seat width options	35, 40, 45, 50 cm	Arm support	Detachable (one-hand operation)
Total weight	29 kg	Seat cushion	Netti Uno pressure redistribution cushion
Max. load	135 kg	Covers	Hygienic Easy Care
Total width	53 - 68 cm dependent on seat width. Seat width +18 cm	Push handles	Height-adjustable push handles
Seat depths	Adjustable by 7.5 cm from 42.5 - 50 cm	Wheels	24" Self propelled main wheels with 7" front castors
Seat heights	46.5/50 cm adjustable	Leg supports	Angle adjustable
Back height	50 cm	Anti-tip	Upwards swingable
Optional back extension	Yes	Head support	Head support A standard. (Height, depth and angle adjustable)
Seat tilt	-5° to +20°	Turning radius	810 mm
Back tilt	90° to 135°	Brakes	Folding user brakes
Back unit	One-piece Velcro back	Accessories	High level of customisation available through addition of accessories
Back cushion	Netti Uno Back	Car Fastening Adapter	

For users that require a specific accessory, option or adaption to configure the Netti 4U CE Plus to their individual requirements, a large range of accessories are available which allows for the chair to be easily customised. This means the chair can be adapted over time to a user's changing needs.

### Available Accessories:

Brake Handle Extension	310 mm	N-86099	Arm Support Pad (Polyurethane)	Large (42 x 8 cm)	N-29137
Smart Back Cushion (with Easy Care Cover)	35x50 cm	N-358304	Hemiplegic Cushion	Right	N-20757
	40x50 cm	N-358306		Left	N-27583
	45x50 cm	N-358308	Arm Support Comfort Pad	Right	N-80917
	50x50 cm	N-358309		Left	N-80918
Stabil Back Cushion (with Easy Care Cover)	35x50 cm	N-368304	Transport Wheel Kit with Drum Brake (pair)	16"	N-29316
	40x50 cm	N-368306		Oval	N-28383
	45x50 cm	N-368308		Square	N-28384
	50x50 cm	N-368309		Head Support (with Easy Care Cover)	Comfort
35x50 cm	N-378301	Adjustable Side Supports	N-90490		
Super Stabil Back Cushion (with Easy Care Cover)	40x50 cm	N-378303	Reinforced Tray	Chin Support	N-91476
	45x50 cm	N-378305		35-43 cm	N-15100
	50x50 cm	N-378306		43-50 cm	N-15101
Pressure Relief Seat Cushion (with Easy Care Cover)	35x50 cm	N-188305	Swingable Tray	35-43 cm	N-15140
	40x50 cm	N-188311		43-50 cm	N-15141
	45x50 cm	N-188317			
Abduction Block	50x50 cm	N-188320			
	Small	N-86096			
	Medium	N-86097			
	Large	N-86098			

\*Ask your Sales Representative about the full range of available accessories for the Netti 4U CE Plus

## A plus benefit for individuality and user friendliness

The Netti 4U CE Plus combines economy and individuality at the highest level. Due to the frame design and the extensive range of accessories, this wheelchair meets various requirements of the user.

A plus benefit for individuality.







1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au