

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

**PRIDE**<sup>®</sup>  
MOBILITY

**GO GO**<sup>®</sup>  
CARBON



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

## Exceptional Portability

The Go Go® Carbon is the perfect blend of strength and portability. Built with a lightweight carbon fibre frame, this 4-wheel scooter offers exceptional durability while weighing just 16kg without its (lithium) battery. It folds down in seconds for convenient transport and storage, and can be towed like luggage through airports and other environments. Powered by an NFC reader and digital display screen the Go Go® Carbon is designed for a lifestyle on the go.

## Specifications

Width	440mm
Length	1023mm
Seat Width	356mm
Seat Depth	275mm
Seat to Deck Height	430mm
Seat to Ground Height	570mm
Ground Clearance	38mm at Motor
Turning Radius	1346mm
Tyres   Front	177.8mm   Solid
Tyres   Rear	203.2mm   Solid
Range per Charge	<16km
Maximum Speed <sup>1</sup>	6kph
Battery Charger	Off-board   2A Lithium ion
Battery Requirements	24V 12 Ah lithium   288Wh
Battery Weight	1.8kg
Weight without Battery	16kg
Heaviest Piece   does not disassemble	16kg
Colour Options	Black   Yellow
SWL	136kg
Warranty	12 months



Conveniently Folds for Easy Portability and Storage



Vibrant Colour Display on Console with NFC Ignition Security



Optional Extra Armrests

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.



Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go Go® Carbon and its lithium-ion battery are permitted.





1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au