

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



*Wonder*TM
ROLLATOR
HD

World's Most Portable Bariatric Rollator!



1300 000 030 | rehabhire.com.au

The World's Most Portable Bariatric Rollator!

The world's most portable heavy-duty bariatric rollator provides reliable support for everyday use at home. It is built for stability and durability on uneven outdoor surfaces, and its strong yet compact design folds easily for travel.

Features and Benefits

- Folds 3x smaller than the average bariatric rollator
- Extra-wide seat and backrest
- Supports larger user weights
- Adjustable handle height
- Large wheels and hand brakes for better mobility
- Basket and organiser pouch for storing personal items
- No assembly required

Specifications

Width	635mm
Depth	685.8mm
Height	787.4mm - 971.55mm
Seat Width	533.4mm
Seat Depth	254mm
Seat Height	558.8mm
Seat Width	533.4mm
Folded Width	406.4mm
Folded Depth	254mm
Wheel Size	203.2mm
Weight	13.5kg
Basket SWL	10kg
Pouch SWL	5kg
SWL	250kg
Warranty	12 months limited



Black Walnut



Regal Rose



Supports up to 250kg



Adjustable handle height



Latch to prevent opening when folded



Locking hand brakes





1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au