

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

PRIDE[®]
MOBILITY

Jazzy[®]
carbon HD



1300 000 030 | rehabhire.com.au



Heavy Duty. Lightweight.

Weighing just 14.5 kg and boasting a 181kg weight capacity, the Jazzy Carbon HD is a game-changing heavy-duty travel power wheelchair. It folds and locks into position in mere seconds and will keep you comfortable with a padded seat back and extra-wide foot platform.

Features and Benefits

- Adjustable lap belt
- Padded back storage pocket
- Durable, lightweight carbon fibre frame
- Compact base folds and unfolds for easy transport
- Height adjustable foot platform
- Latch for locking unit open or closed



Full Front and Rear Suspension



Quick Disconnect Joystick
with LCD Display



Puddle Light



Specifications

Width ²	670mm
Length ²	935mm
Seat Width	510mm
Seat Depth	510mm
Seat to Ground Height	500mm Fixed
Ground Clearance	25mm at Motor
Turning Radius	925mm
Front Castors	165.1mm Solid
Tyres Rear	269.24mm Carbon Fibre
Range per Charge ¹	<20km
Maximum Speed ¹	6kph
Battery Charger	Off-board 3.5A Lithium-ion
Battery Requirements	20Ah, 24V Lithium 436.8Wh
Battery Weight	3.4kg
Heaviest Piece	14.6kg
Colour	Black Cherry Carbon
SWL ³	181kg
Compliance	ARTG 342775 AS NZ 3695.2:2019
Warranty	6 12 24 months limited



1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Width of base (not including seat).
3. Maximum weight capacity is reduced to 181 kg on 510mm wide seat configuration.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

Notes



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au