

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

PRIDE[®]
MOBILITY

Jazzy[®]
ULTRA LIGHT



1300 000 030 | rehabhire.com.au



Your Greatest Travel Companion

Experience the ultimate in lightweight performance and strength with the Jazzy® Carbon Ultra-Lite powerchair. Designed for active lifestyles, it features a durable carbon fibre frame that folds down in a few quick motions and is powered by a reliable lithium-ion battery. The Jazzy® Carbon Ultra-Lite provides effortless mobility and control, combining convenience with lightweight durability.

Features and Benefits

- Two motor, 24v, 170W Rear wheel drive
- Bright LED lighting for visibility in all environments
- Integrated foot platform
- Front castor forks for enhanced manoeuvrability
- Lightweight, durable carbon fibre frame
- Long-lasting, lithium-ion battery
- Compact one lever fold function

Features and Accessories



Long-lasting, Lithium-ion Battery



Rear LED Lighting



Compact One Lever Fold Function



Travel Bag



Mesh Cup Holder



Armrest Padding



Attendant Mount



Headrest

Specifications

Width ²	620mm
Length ²	950mm
Seat Width	450mm
Seat Depth	450mm
Seat to Ground Height	480mm Fixed
Ground Clearance	55mm at Motor
Turning Radius	900mm
Front Castors	203.2mm Solid
Tyres Rear	304.8mm Carbon Fibre
Range per Charge ¹	<20km
Maximum Speed ¹	6kph
Battery Charger	Off-board 3.5A Lithium-ion
Battery Requirements ³	10Ah, 25.2V Lithium 252Wh 20Ah, 25.2V Lithium 504Wh
Battery Weight	10AH: 2.3kg (standard) 20AH: 3.2kg (optional)
Weight without Battery	16kg
Heaviest Piece ⁴	15.3kg
Colour	Carbon
SWL	136kg
Compliance	ARTG 342775 AS NZ 3695.2:2019
Warranty	6 12 24 months limited



1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Width of base (not including seat).
3. Standard 10.4 Ah battery is airline compliant. Optional 20 Ah battery is not safe for air travel.
4. Seat weight included. Battery weight not included.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth adjustable.



 Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go Go® Carbon and its lithium-ion battery are permitted.

Notes



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au